


1527



-lektivid	337.17	1.115 · 10 ⁻² ml	3.76 g
dioxanon	130.15	2.230 · 10 ⁻² ml (2 eq)	2.9 g
R-prolin	115.13	3.34 · 10 ⁻³ ml (30 ml / 1)	385 mg
DMF	22.5 ml		

Jon. 17⁰⁰ Poate 1 dena → videti se proizvod, di ima dosta ASD - shodato 1 ml dioksanona i 100 mg proline jso. Ostajemo de la selje mezz. Poate 2 dena - jso 1 ml dioksanona - poate 3 dena jso 1 ml dioksanona - petele nijstun - gotava. Obzaten - razobica chun i vodom, vodeni dno elupolovan 2x (50+200 ml sh0) i spogim se prvim sh0 (do 20 ml) → spogim ozyen shi spogim ipem 1x 40 (100 ml), bime-on on rano i upremo. Dry flash 185 g diol, elut PE/ET = 85/15, izasla od (shuro ~ 8g) (trahoj do 100 ml)

29 - 51 eq tuc, izluceno 3.15g (61%).

